MutaCHIP[®] TOXO

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MutaCHIP[®] TOXO

DNA-Macroarray Test Kit for the Analysis of Genetic Variations
Associated with Detoxification

For in vitro diagnostics only

KF391011

10

KF390011

20

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany
www.immundiagnostik.com
info@immundiagnostik.com
 Tel.: +49 (0)6251/ 701900
 Fax: +49 (0)6251/ 849430

Content

1 Application	3
2 Introduction	4
3 Concept of the Assay	4
4 Components	5
5 Materials	5
6 Storage and Shelf life	6
7 Working Conditions	6
8 Considerations and Precautions	6
9 Sample Collection	6
10 Test Procedure	7
1 PCR Product Generation.....	7
2 PCR Protocol.....	8
3 MutaCHIP Protocol.....	9
11 Evaluation and Interpretation of Results	10
12 Trouble shooting	14

1 Application

The MutaCHIP[®] TOXO kit is a biomolecular test for the detection of mutations in genes responsible for the detoxification of the body. The genetic variations examined are associated with the biotransformation of xenobiotics and the antioxidative system. Genomic DNA is used for the analysis by means of a DNA macroarray. The following single nucleotide polymorphisms (SNP) and mutations are analysed by the MutaCHIP[®] TOXO:

Character	Gene	Variation	rs number
Phase I Biotransformation	CYP1A1	*2A (3798T>C)	rs4646903
	CYP1A2	*1C (-3860G>A)	rs2069514
		*1F (-163A>C)	rs762551
	CYP2C9	*2 (430C>T)	rs1799853
		*3 (1075A>C)	rs1057910
	CYP2C19	*2 (19154G>A)	rs4244285
		*3 (17948G>A)	rs4986893
		*17 (-806C>T)	rs12248560
	CYP3A5	*2 (27289C>A)	rs28365083
		*3 (6986G>A)	rs776746
Phase II Biotransformation	GSTM1	14kb Deletion	
	GSTT1	50kb Deletion	
	GSTP1	A>G (Ile105Val)	rs1695
	NAT2	481C>T	rs1799929
		590G>A	rs1799930
		857G>A	rs1799931
	COMT	G>A	rs4680
Phase III Biotransformation	MDR1	3435C>T	rs1045642
Antioxidative System	SOD2	C>T	rs4880
	CAT	-262C>T	rs1001179
	CYBA	242C>T	rs4673
	NOS3	-786C>T	rs2070744
Efficacy and Side-Effects of Warfarin	VKORC1	*2 (-1639G>A)	rs9923231

2 Introduction

The human body is continuously confronted with foreign and to some extent toxic substances. These xenobiotics include drugs, stimulants and chemicals. A great many of these foreign substances cannot be readily excreted as they are lipophilic, i.e. poorly or not at all water-soluble. Biotransformation prevents accumulation of xenobiotics by transforming lipophilic compounds into more hydrophilic ones that can be excreted from the body. Additionally biotransformational processes are involved in converting pro-drugs (inactive medication precursors) into their active forms.

There are three phases of biotransformation. During phase I functional groups are attached, either via oxidation, reduction or hydrolysis. Whilst phase II the intermediate products from phase I are coupled to endogenous, mostly highly hydrophilic functional groups (conjugation reaction). The last phase, phase III, encompasses export of now water-soluble substances from the cells. If this process is hampered by enzymatic deficiencies, e.g. caused by polymorphisms and other genetic variations, severe complications can develop. Foreign substances accumulate in the body potentially leading to serious side-effects or diseases, as for example cancer. Consequently, also pro-drugs are not activated and are not fully functional and efficient.

Reactive oxygen species (ROS) are damaging to all kinds of biomolecules. During biotransformation, as well as by drugs and environmental toxins, ROS can be produced inside the body. The antioxidative system is engineered to eliminate ROS and protect cells before they are damaged. Among others, cancer and neurodegenerative diseases can be the consequence of impaired ROS clearing, e.g. by dysfunctional enzymes.

A number of polymorphisms and gene variations with relation to detoxification have been identified so far (see section 1). Genotyping detects these mutations and identifies such persons who are at a higher risk for developing the above mentioned conditions.

References:

- Coleman M. D. Human Drug Metabolism: An Introduction. 2. Edition, 2010. Hoboken, NJ, USA: Wiley-Blackwell.
- Ernst B., Vögtli A. Moderne Pharmakokinetik: Transport durch Membranen. 1. Edition, 2010. Weinheim, Germany: Wiley-VCH Verlag GmbH & Co KGaA.
- Karlson P., Doenecke D., Koolman J., Fuchs G., Gerok W. Karlsons Biochemie und Pathobiochemie. 5. Edition, 2005. Stuttgart, Germany: Thieme Verlagsgruppe.

3 Concept of the Assay

To analyse the mutations a PCR is performed, which amplifies the relevant genes using freshly extracted genomic DNA as template. The amplification product is transferred to the MutaCHIP, where it binds to the immobilized probes. During the washing step not perfectly bound fragments are removed. Subsequently, the enzyme is added, which binds to the remaining probe-fragment complexes. After adding the substrate the precipitation reaction occurs only at those spots at which DNA is still attached. The colored precipitation pattern is then detected by the image reader and analyzed as well as interpreted by the associated software. The analysis via precipitation reaction allows a distinct genotyping of all alleles spotted on the MutaCHIP.

4 Components

Each test kit contains the following components for 10 / 20 MutaCHIP assays and the instruction manual:

Description	Size of Reaction Tube or Flask	Number
PCR Mix A	2 mL Reaction Tube (green lid)	1
PCR Mix B	2 mL Reaction Tube (yellow lid)	1
PCR Mix C	2 mL Reaction Tube (red lid)	1
PCR Mix D	2 mL Reaction Tube (blue lid)	1
PCR Mix E	2 mL Reaction Tube (white lid)	1
Polymerase	2 mL Reaction Tube (purple lid)	1
Hybridisation Buffer	15 mL / 30 mL Plastic Flask	1
Washing Buffer	30 mL / 60 mL Plastic Flask	1
Conjugation Mix	2 mL Reaction Tube (black lid)	1 / 2
Substrate	2 mL Reaction Tube (brown, blue lid)	1 / 2
DNA Ref	2 mL Reaction Tube (orange lid)	1
MutaCHIPs	1.5 mL Reaction Tube	10 / 20
Instruction manual	/	1

5 Materials

Required materials that can be ordered separately from PharmGenomics:

- Macroarray Analysing Tools (MAAT):
 - Notebook + TOXO Genotyping Software
 - MutaCHIP image reader
 - Thermocycler (Pqclab Primus 25 advanced)
 - Thermoshaker with cooling function (BIOR Mixing Block MB-102)

Required Materials – not provided

- pipettes:
 - 0.1 - 2.5 µL
 - 0.5 - 10 µL
 - 10 - 200 µL
 - 100 - 500 µL
- 0.2 mL PCR tubes (sterile)

6 Storage and Shelf life

- The Polymerase has to be stored at -20 °C
- All other components are stored at 2-8 °C.
- The substrate has to be strictly protected against light exposure.
- The bags contain each 5x MutaCHIPs with open lids. After unsealing a bag, the lids of the remaining MutaCHIPs have to remain open as the protection gas escapes.
- The MutaCHIPs can be stored in the loosely (not airtight) closed (do not use tape) light-protection bag up to several weeks and room temperature (RT). For storage choose a dry and dark place.
- To avoid even minimal loss of performance, we recommend using up one bag of MutaCHIPs within two weeks.
- The MutaCHIPs have to be protected against direct exposure to sunlight and dust. and dust.

7 Working Conditions

- Never centrifuge the MutaCHIPs.
- Do not touch the surface of the MutaCHIPs with a pipette.
- Use only substances that are mentioned in the protocol.

8 Considerations and Precautions

The guidelines and principles for working in a biomolecular laboratory have to be followed.

- The test is suitable for genomic DNA **freshly** extracted from whole EDTA blood as starting material. Only then optimal results are guaranteed.
- Do not mix reagents from different lots.
- The MutaCHIPs are:
 - for single-use only
 - only for *in vitro* diagnostics
- Do not let the MutaCHIP dry out while performing the analysis!
- The MutaCHIP is designed for the use with the MAAT and the software provided by PharmGenomics.
- Perform all steps in a timely manner.
- Keep all stock solutions cooled while working.
- Always open the MutaCHIP with both hands. Do not apply pressure to the tube.

9 Sample Collection

The template for PCR amplification is genomic DNA from EDTA whole blood. The DNA concentration should be between 15 and 30 ng/μL. The minimum DNA purity (A260/A280 ratio) should be higher than OD 260/280 1.8.

For the assay high molecular (freshly extracted) DNA has to be used.

10 Test Procedure

10.1 PCR Product Generation

- Before starting, make sure that all working surfaces and utilized devices are decontaminated.
- Prepare all components of the kit required for the procedure of the PCR
- Subsequently, perform the PCR according to the protocol listed below.
- For the amplification of target DNA, five PCR reactions per sample are required - one for each mastermix! Those are pipetted as described in the following tables.

PCR Mix A:

Component	Volume per 25 µL- reaction mixture
DNA (min. 60 ng - max. 120 ng)	4.0 µL
PCR Mix A (green lid)	20.8 µL
Polymerase	0.2 µL

PCR Mix B:

Component	Volume per 25 µL- reaction mixture
DNA (min. 60 ng - max. 120 ng)	4.0 µL
PCR Mix B (yellow lid)	20.8 µL
Polymerase	0.2 µL

PCR Mix C:

Component	Volume per 25 µL- reaction mixture
DNA (min. 60 ng - max. 120 ng)	4.0 µL
PCR Mix C (red lid)	20.8 µL
Polymerase	0.2 µL

PCR Mix D:

Component	Volume per 25 µL- reaction mixture
DNA (min. 60 ng - max. 120 ng)	4.0 µL
PCR Mix C (blue lid)	20.8 µL
Polymerase	0.2 µL

PCR Mix E:

Component	Volume per 25 µL- reaction mixture
DNA (min. 60 ng - max. 120 ng)	4.0 µL
PCR Mix C (white lid)	20.8 µL
Polymerase	0.2 µL

Place the samples in the Thermocycler and continue with chapter 10.2.

10.2 PCR Protocol

The PCR protocol has to be established anew for each laboratory, as different thermocyclers have different heating rates. This establishment can be omitted when using the recommended thermocycler (Peqlab Primus 25 advanced). The following standard protocol can be used to start the establishing process:

	Temperature [°C]	Time [min]	
Start	94	2	1x
Denaturation	94	0.5	10 x
Primer Hybridisation	72	0.5	
Elongation	72	1.5	
Denaturation	94	0.5	25 x
Primer Hybridisation	60	0.5	
Elongation	72	1.5	
Final Elongation	72	5	1x
Store at	8	Forever	-

10.3 MutaCHIP Protocol

A) Preparation of the Hybridisation Buffer

If the Hybridisation Buffer is turbid or a precipitate can be seen, it has to be heated in a microwave for several seconds (240 W) or in a water bath. Gently shake the Hybridisation Buffer until it gets clear, to homogenise it again. Before usage, it has to be cooled to room temperature (RT).

B) Preparation of DNA samples

- Preheat the thermoshaker to **51 °C**
- Add **5 µL** PCR product A, **5 µL** PCR product B, **5 µL** PCR product C, **5 µL** PCR product D und **5 µL** PCR product E in a separate tube and mix by brief vortexing
- Denature the mixture for **2 min** at **95 °C** in a thermocycler
- After denaturation, add **100 µL** Hybridisation Buffer
- Pipette the denatured mixture into the MutaCHIP without touching the bottom of the tube

C) Hybridisation

- Perform the hybridisation at **51 °C**, **550rpm** for **30 min**.

D) Washing steps after Hybridisation

Caution: Pay attention to the correct adjustment of volume!

- Completely remove the Hybridisation Buffer from step C)
- Add carefully **500 µL** of Washing Buffer onto the MutaCHIP
- Wash the MutaCHIP at **550 rpm** and **51 °C** for **5 min**

E) Conjugation step

Set the thermoshaker to **21 °C**.

Note: During the cooling period the MutaCHIP has to be removed from the thermoshaker! The Washing Buffer has to remain on the MutaCHIP until the target temperature is reached!

Caution: Pay attention to the correct adjustment of volume!

- Completely remove the Washing Buffer from step D)
- Add **100 µL** of Conjugation Mix to the MutaCHIP
- Incubate the MutaCHIP at **550 rpm** and **21 °C** for **15 min**

F) Washing step after conjugation step

Caution: Pay attention to the correct adjustment of volume!

- Completely remove the Conjugation Mix from step E)
- Add of **500 µL** Washing Buffer to the MutaCHIP
- Wash the MutaCHIP at **550 rpm** and **21 °C** for **5 min**

G) Precipitation

Caution: Pay attention to the correct adjustment of volume! **Do not shake the MutaCHIP during the precipitation!**

- Completely remove the Washing Buffer from step F)
- Add **100 µL** of Substrate to the MutaCHIP and incubate for **5 min** at **21 °C**. To do so, put the MutaCHIP into the thermoshaker (**Do not activate shaking function**)
- Thereafter, remove the Substrate completely (pipette) and immediately add **500 µL** of Washing Buffer
- Place the MutaCHIP into the image reader (pay attention to the correct orientation), take a picture and start with the analysis/evaluation of the results (see following chapter)


11 Evaluation and Interpretation of Results

The evaluation is performed using the MAAT and the TOXO Genotyping Software. The results are compiled with help of the software into a report. For the evaluation of the chip, follow the short instructions below. For further and detailed information, please refer to the TOXO Genotyping Software handbook.


Step 1: Create a new project

Click on the button *New experiment*. Assign an arbitrary name for the experiment and subsequently save it by clicking on the button *save*.


Step 2: Start the analysis process


Click on the *Start* button to initiate data analysis.

Step 3: Quality check of the MutaCHIP

To ensure a correct analysis result, the picture quality of the MutaCHIP has to be checked. Particles of dust on the bottom side of the MutaCHIP can affect the analysis. These can be removed by cleaning with a soft and wet tissue. Follow the instructions given in the screen shown here. Press *Continue analysis* if the quality of the picture is comparable to the figure above.

Step 4: Genotyping results

After the complete data analysis, the results can be accessed in the analysis module / genotyping module or in the diagnostic report (example picture). The used icons are explained in the following table.

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
	The patient carries on both alleles the healthy variant of the investigated genetic variants.
	The patient carries on one allele the healthy and on the other allele the mutated genetic variant.
	The patient carries on both alleles the mutated genetic variation in homozygous form.
	The signal values of the probes for this genetic variation are too weak for a valid result. This could be caused by other sequence variations in close proximity to the investigated variant. The remaining signals of the assay are not influenced.
	Both probes for this genetic variant show an invalid signal. This might be caused by impurities or dust particles on the MutaCHIP surface. The remaining signals of the assay are not influenced.

Step 5: Diagnostic report

To evaluate the results open the diagnostic Report Main. Therefore click on the button **Report**.

Additionally a .pdf document can be created or the report can be directly printed.

12 Trouble shooting

Problem	Solution
Software Message: Warning! The signals of the biotin reference markers are too low. This could be a sign of a failed or missed conjugation step. Alternatively the enzyme could be degraded. The system will be stopped.	Check enzyme and/or substrate. Repeat the assay with new enzyme/ substrate
Software Message: Warning! The signals of the DNA probes indicate a failed amplification.	Check the PCR products on a 2% agarose gel. If the gel is empty, please perform the amplification again and check on a gel before performing a chip.
No bands on agarose gel	Perform the DNA extraction again. If the problem persists, check the PCR protocol according to chapter 10.2 or contact PharmGenomics customer service.
Software Message: No error description ErrorCode -3011	The reader is not plugged in properly: Keep "Esc" pressed for 3 seconds and plug in the reader in the right manner – please refer to the user manual of the PGDx System
Poor image quality	Clean the bottom side of the MutaCHIP using a cotton swab or a cloth wet with disinfectant Repeat the photograph (it can be necessary to repeat this procedure several times)
Blurred picture	Carefully clean the camera with a cotton swab